

# Anti Choking Trainer



## [Product Description]

Model: WAP-ACT-01

Accessories: 1pcs Vest, 1pcs user manual

1pcs bag, 1pcs foreign body obstruction

Material: Diving Fabric


N. W.: 1.2kgs

Usage: For training use.

Warranty: 1 year.

The Anti Choking Trainer is worn by students in CPR class to learn the Abdominal Thrust Maneuver. When performed correctly, a foam plug shoots into the air, making instruction realistic and leaving students confident. This is the easiest and most comfortable device to teach choking skills on the market. It is truly a must for any CPR instructor!

## How To Use ?



- A.** Put the foreign body into the throat of the inlet. Stand or kneel behind the patient, wrap your hands around the patient's waist, and make a fist with one hand.
- B.** Put the held thumb against the patient's abdomen, on the midline of the abdomen above the navel and below the breastbone.
- C.** Hold the holding hand with the other hand and quickly hit the patient's abdomen upwards. Repeated rapid shocks guide the expulsion of foreign objects from the airway.
- D.** Use the round back pad for slap training.