



A501

ELECTRIC-HYDRAULIC OPERATING TABLE

Registration code: SKT-B



SAIKANG/Medical

Healthier
your life

A501

FEATURES AND HIGHLIGHTS



Technical parameters

▪ External size (LxWxH):	2125x550x(600-950)mm
▪ Reverse trendelenburg:	0-15°
▪ Trendelenburg:	0-20°
▪ Lateral tilt left/right:	0-15°
▪ Head plate up/down:	40°/90°
▪ Back plate up/down:	85°/40°
▪ Leg plate up/down:	15°/90°
▪ Kidney bridge elevation:	0-800mm
▪ Horizontal sliding:	0-360mm
▪ Horizontal rotation:	180°
▪ Flex:	95°
▪ Reflex:	220°

Technical configuration

☒ Bed frame	1pc
☒ Mattress	1set
☒ Anesthesia screen frame	1pc
☒ Shoulder holder with cushion	2pcs
☒ Arm holder with cushion	2pcs
☒ Leg support with cushion	2pcs
☒ Foot support with cushion	2pcs
☒ Hand controller	1pc
☒ Kidney bridge handle	1pc

Material

The base column and side holder are all made of high quality 304 stainless steel, which is anti-corrosive and can be easily cleaned.

**Dual control system**

All functions controlled by both hand controller or control panel on table base column. With action LCD display on hand controller, easy to operate. With control panel on table base column, in case of emergency hand controller failure.

**Special functions**


With table top 180° rotation function.

**Special functions**


With one-key reset, one-key flexion and re-flexion function.

**Brake**




Electric Control system for the brake, it is easily for moving operating table.

**Kidney bridge**

Built-In Kidney Bridge, More Convenient And Nice Appearance.

**Mattress**

PU cover, memory mattress, easy cleaning, thickness of mattress: 5cm.

**Anesthesia screen frame****Leg support with cushion
Foot support with cushion****Arm holder with cushion****Shoulder holder with cushion****Body holder with cushion**

Professional functions

Lowest position of table top

600mm

Highest position of table top

950mm

Back plate up/down

85°/40°

Kidney bridge elevation

0-80mm

Reversed trendelenburg

0°-15°

Trendelenburg

0°-20°

Right lateral tilt

0°-15°

Left lateral tilt

0°-15°

Head plate upward

0°-40°

Head plate downward

0°-90°

Leg section evagination

0°-90°

Leg plate down

0°-90°



Design



Research



Manufacturing



Maintenance



Service



South America Branch

Saikang Medical Technology, Inc
W3400 Corrientes, Argentina
TEL: +55 (11) 2758-0438
Rua Universal, Hollywood,
São Bernardo do Campo,
Brazil
TEL: +55 (11) 2758-0438

East Africa Branch

Saikang Medical Technology Ltd
Mezzanine Floor, Taj Tower,
upper hill, nairobi, Kenya
TEL: +254714084806

South East Asia Branch

Saikang Medical Technology Ltd
Sala Thammasop Subdistrict
Thawi
Watthana District
Bangkok Thailand
TEL: +66 642741046

UK Branch

Uk Healthward International Limited
DEPT 111,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +31 (0) 202111106

Russia Branch

LLC 'SAIKANG'
454004, Russian Federation,
Chelyabinsk,
Universitetskaya Naberezhnaya,
Russia
TEL: +7 (919)300-00-06

Mid East Branch

Saikang Medical Supplies (Mid East) LLC
Musaffah
P.O Box 53200, Abu Dhabi,
UAE
Tel: +971 585710652

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎ +86 18021231901 · 🕒 24/7
✉ export@saikangmedical.com
🌐 www.saikangmedical.com



Company video



General catalogue