

BESDATA® BD-M7DP Disposable Video Laryngoscope with 7 Inch Monitor

1

Single use blades, reducing the risk of cross-infection, reprocessing costs, time and storage.

2

Anti-Fog Design

3

Seperate 7inch LCD Monitor, to keep doctor away from patients secretions.

4

More Available Space with 12mm BESDATA blades height

5

Blade work with both BD-M3 portable monitor and BD-M7 7inch Monitor



Model Number: BD-M7DP

Order number	BD-M7	DP-MiI0	DP-MAC1	DP-MAC2	DP-MAC3	DP-MAC4
Description	7in LCD Monitor	MiI0	MAC1	MAC2	MAC3	MAC4
Application	LCD Monitor	Neonate	Pediatric	Small Adult	Adult	Large Adult
Blade Length(cm)	/	8.31	10.15	12.11	16.08	17.5

Technical Data

Parts	Description	Technical Index
Monitor	Display	3 inch LCD
	Resolution Ratio	640*480(RGB)
	Viewing angle	>60°
	Power	<2W
	Display front and back rotation angle	0°~130°
	Monitor Thickness	12mm
	Weight	<90g
Camera	Illuminance	>150 LUX
	Resolution Ratio	1280*720 px
Image/Video function	Image/ Video function	Yes
	Output	USB output, Easy to storage and establish image /video
	Memory type	Micro SD card 8GB
	Video/Photo Play back	Yes
Battery	Battery Type	Rechargeable lithium battery
	Capacity	730mAh
	Cycle life of battery	>300 times
	Battery working time	>75min
	Charging time	<90min
Blades	Weight	< 80g