

BESDATA® BD-DP Disposable Video Laryngoscope



1

Single use blades, reducing the risk of cross-infection, reprocessing costs, time and storage.

2

Anti-Fog Design

3

The ergonomic design, easy to intubation, easy to hold.

4

More Available Space with 12mm BESDATA blades height

5

Blade work with both BD-M3 portable monitor and BD-M7 7inch Monitor



Model Number: BD-DP

Order Number	BD-M3	DP-Mil0	DP-MAC1	DP-MAC2	DP-MAC3	DP-MAC4	BD-DP-H1	BD-DP-H2
Description	3inch LCD Monitor	Miller 0	MAC1	MAC2	MAC3	MAC4	Neonate Handle For Miller 0 blade	Adult Handle for MAC1, MAC2, MAC3, MAC4 blades
Application	LCD Monitor	Neonate	Pediatric	Small Adult	Adult	Large Adult		
Blade Length	/	8.31	10.15	12.11	13.69	16.08		



Technical Data Sheet

Parts	Description	Technical Index
Monitor	Display	3 inch LCD
	Resolution Ratio	640*480(RGB)
	Viewing angle	>60°
	Power	<2W
	angle	0°~130°
	Display up and down rotation angle	0-270°
Camera	Illuminance	>150 LUX
	Resolution Ratio	1280*720 px
Image/Video function	Image/ Video function	Yes
	Output	/video
	Memory type	Micro SD card 8GB
	Video/Photo Play back	Yes
Batterty	Battery Type	Rechargeable lithium battery
	Capacity	1350mAh
	Cycle life of battery	>300 times
	Battery working time	>140min
	Charging time	<150min