

SP5 SYRINGE PUMP

- Compact & Portable, easy to use
- LED and LCD dual display
- Dual CPU system, safe & more reliable
- Manual bolus and automatic bolus
- Flow rate, time/volume, dosage/weight modes
- All standard syringes applicable including 5, 10, 20, 30, 50/60ml
- Unique direct pressure sensor, dynamic pressure value and picture display, 3-level occlusion sensitivities setting
- Particular built-in sliding potentiometer ensures automatic calibration and accuracy
- Automatically release impact dose after occlusion (anti-bolus)
- Vertical & Horizontal rotatable clamp
- 2000 history records and last configuration memory



Global Markets



Amoul Ambulanc (Shenzhen) Tech. Co., Ltd.



Building 5C, Skyworth Innovation Industry Park, Tangtou 1st Road,
Baoan District, Shenzhen, P.R.China
+86 755-26072286
www.amoulmed.com
info@amoulmed.com

Amoulmed ambulshenzhen.Tech Amoulmed Amoulmed Amoul

© Copyright 2022 Ambulanc(Shenzhen) Technology Co., Ltd
SZAB IPC/E20200507-V1.0 Printed in China

INFUSION SYSTEMS

Safe Efficient Reliable



IP5 Infusion Pump

Comprehensive Infusion Modes to Ensure Patient Safety

- Compact & Portable
- Ergonomic intuitive interface, easy to operate
- Allow to change the flow rate while operating
- 3.2-inch high-definition color display
- Dual CPU system, safer & more reliable
- Manual bolus and automatic bolus
- Flow rate, time/volume, dosage/weight modes
- Different infusion mode in a particular stage
- Drip sensor and fluid heater (Optional)
- All brands infusion sets applicable
- Dynamic pressure display, 3-level occlusion sensitivities setting (the value of each level is programmable)
- Automatically release impact dose after occlusion (anti-bolus)
- Vertical & Horizontal rotatable clamp
- Drug library, wireless monitoring upgradable



SP6 Double Channel Syringe Pump

The Safest & Most Efficient Dual CPU System

- Operation interfaces separated without interfering with each other
- Dual CPU system, safe & more reliable
- Manual bolus and automatic bolus
- Allow to change flow rate while running
- Flow rate, time/volume, dosage/weight modes
- All standard syringes applicable including 5, 10, 20, 30, 50/60ml
- 3-level occlusion sensitivities setting
- Particular built-in sliding potentiometer ensure automatic calibration and accuracy
- Automatically release impact dose after occlusion (anti-bolus)
- Vertical & Horizontal rotatable clamp
- 2000 history records and last configuration memory
- High-definition color display
- Drug library, wireless monitoring upgradable

