

MAICO MA1 HEARING SCREENER WITH HEADSET



Product Code	33600
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	2233027
CND	Z121401
EAN/UPC	4260176120858

Description

MAICO MA1 HEARING SCREENER WITH HEADSET

Hearing screener with super-lightweight design, small enough to fit in the palm of the hand. Ideal for fast screening purposes. The device is capable of delivering a pure tone signal to one ear at a time. The stimulus range is 15 to 50 dB HL in level steps of 5 dB. Test frequencies range from 500 to 4000 Hz. Automatic ON/OFF function. It is supplied with headset DD45, cradle and 2AA batteries. Multilingual manual: GB, FR, IT.

Technical Specifications

Test signal: pure tone (continuous)
Test frequencies: 500, 1000, 2000, 4000 Hz
Level steps: 5 dB
Hearing level range: 15dB HL to 50 dB HL
Accuracy - distortion:± 4 dB - less than 3%
Size - weight: 6.4x15.3xh 2.3 cm - 115 g
Power supply: 2AA batteries (included)
Standards: IEC 60601-1:2012 Type B, IEC60601-1-2:2014
Classification: Class IIa (EU 2017/745) Medical Device