

Lepu Medical Technology (Beijing) Co., Ltd was established in 1999. It is specialized in developing, manufacturing and marketing high-tech medical devices and equipment. Today, Lepu Medical has grown into a global leading medical technology company in the fields of cardiovascular interventions, structural heart diseases, cardiac rhythm management, anesthesia and critical care, in vitro diagnostics and general surgery. Our vision is to develop Lepu Medical into one of the largest supplier on the platforms of medical device, medicine, healthcare products, healthcare service and new medical mode.



Lepu Medical Technology (Beijing) Co., Ltd

Add: Building 7, No. 37 Chaoqian Rd., Changping District, Beijing, China 102200
Tel: + 86 755 2643 3514 Fax: + 86 755 2643 0930 E-mail: info@creative-sz.com
<https://en.lepumedical.com>



iHolter

Real-time Holter Monitor

Realize AI-ECG + Real-time Monitoring
+ Regional ECG Network



Continuous
Tracking



Risk Warning



Black Box



AI-ECG Analysis



Remote
Monitoring



Professional
Report

Science and technology cherish life



Product Features



- Lightweight and portable, 2.4" color screen
- Support lead detachment alert and low battery alert
- Support 4G/USB data transmission
- Support pacemaker detection
- One key to record when feeling uncomfortable

Models	iHolter-3 with 6 electrodes	iHolter-3 with 10 electrodes	iHolter-12 with 6 electrodes	iHolter-12 with 10 electrodes
Display of Waveforms on the screen	3 waveforms	3 waveforms	6 waveforms	12 waveforms
Display of Waveforms on AI-ECG Tracker	8 waveforms	12 waveforms	8 waveforms	12 waveforms
Working Time	7 days	5 days	7 days	5 days
Battery	Users should alternate between 2 batteries, 1 battery supports 30 hours			

Real-time Monitoring Solution



AI-ECG Tracker (Cloud Version)

- Support real-time monitoring
- Support real-time display of HR
- Provide abnormal ECG event warnings
- Support display of record progress

Traditional Holter Solution

No network required for data export and import and support different deployment versions of AI-ECG Tracker (including stand-alone, local and cloud version).


Record Data → Import Data → Data Analysis → Generate Reports



Advanced AI-ECG Analysis

AI-ECG Tracker is a professional analysis & data management software for holter monitor.

- **Improve accuracy:** comprehensive accuracy rate of 95.2%
- **Reduce analysis time:** analysis within 5 minutes (24h data)
- **Improve clinical classification:** 39 types of abnormal events



Clinical Application

