

SEESHEEN

Cy

 Image

 **CY-E11 Series**
FLEXIBLE VIDEO CYSTOSCOPE

 0197

Tapered Tip

Tapered tip design for easier insertion.



Multi-function Buttons

- ① Zoom in/out
- ② Capture/record
- ③ Freeze/release



Angle Lock Design


Lock or release the angle of the distal end for accurate diagnosis and treatment.



*Patent NO.: ZL202022141661.X

Suction Function

One-button suction for clear view. Unique anti-drop design to prevent the suction valve from falling off.



*Patent NO.: ZL202022141662.4

Specifications

Model	Insertion Tube	Distal End	Working Channel	Working Length	Field of View	Depth of field	Bending (Up/Down)	Rotation (Left/Right)	Suction Function
CY-E11 plus	15.6Fr	14.7/16.2Fr	6.6Fr	380mm	120°	3-100mm	220°/130°	120°/120°	Yes
CY-E11	15.6Fr	14.7/16.2Fr	6.6Fr	380mm	120°	3-100mm	220°/130°	120°/120°	No

Image Processor



10.1"
Touch screen



1280*800
HD Resolution



HDMI+ AV
output



Rechargeable battery
supporting 4 working-hour



Detachable
SD card



Portable
design

With Suction Function



Without Suction Function

