

Toby Bed

02384 The Toby baby bed consists of a metal frame and a perforated metal mattress platform. The sides are made of securit glass.

Frame

- Frame of round steel tubes.
- 4 design castors Ø 75 mm with brakes.

Mattress platform

- Removable perforated steel mattress platform.
- The mattress platform can be adjusted in the reverse-trendelenburg position by means of 2 telescopic adjusts (*).

Sides

- Made of securit glass.
- The sides are firmly attached in the frame.
- Thermally shaped plexiglass cupola on top of the bed.

Colours


We would like to refer to our extensive colour gamut.
Sample card + samples available on request.

(*) Option



Dimensions

- LxWxH | 80 cm x 43,5 cm x 103 cm
- Recommended mattress dimension:
75 cm x 39 cm x 6 cm



Materials and finishing

- Metal: Steel, epoxy-coating or chrome.
- Securit glass.
- Withstands the common cleaning products.

For more information about the **materials** and **construction**: please consult our material lists.

Options and accessories

- Reverse-trendelenburg, adjustable mattress platform.

