

Nino Intensive

08878 A very reliable cot with 4 lowerable sideguards to guarantee an optimal accessibility.

Frame

- The frame is made of round steel tubes.
- Equipped with 4 integrated IV rod cases and setting knobs.
- The bed frame rests on 4 design castors Ø 100 mm with brakes.

Mattress platform types

- **1. Fixed mattress platform:** consists of a metal frame and a solid laminate plate, thickness 6 mm.
- **2. Inclined mattress platform:** consists of a metal frame with a removable solid laminate mattress platform plate. Mattress platform inclination is possible up to $\pm 12^\circ$ by means of 2 telescopic adjusts.
- **3. Inclined 2-section mattress platform:** consists of a metal frame with 2 removable solid laminate mattress platform plates. Inclined headrest up to $\pm 70^\circ$. Both mattress platform sections are equipped with a flexible handle that can only be operated on the outside of the bed.

Sideguards

- The sideguards are made of round steel tubes.
- The distance between the round rails is 60 mm at the most.
- The sideguards can be put in 3 different positions: lowest position, middle high position (33,5 cm) and highest position.
- The lowerable sideguards glide over 2 chromium-plated solid steel guide rails.
- The sideguards can only be unlocked by using your both hands.

The cots are made in accordance with the European norm EN 716 and comply with the Medical Devices Regulation EU 2017/745

Colours

We would like to refer to our extensive colour gamut. Sample card + samples available on request.


Accessories

We would like to refer to our accessory list.



Dimensions

- **Type 55 x 110 cm:**
a x b: 65 x 110 cm
recommended mattress dimension: 55 x 100 x 8 cm
- **Type 70 x 157 cm:**
a x b: 80 x 167 cm
recommended mattress dimension: 70 x 157 x 8 cm



* dimensions for a fixed mattress platform

Materials and finishing

- **Metal:** Steel, epoxy-coating or chrome.
- **Synthetic:**
- Withstands the common cleaning products.

For more information about the **materials and construction:** please consult our material lists.

Options

- Inclined mattress platform.
- Inclined 2-section mattress platform.
- Fire-retardant anti-urine mattress.
- IV rod to install in the cases.
- Double oxygen bottle holder.
- Shelf/monitor shelf.
- Linea twin castors Ø 100 mm.

