

CloudHolter

Ultra-light and reliable recorder

High quality ECG recording

Powerful and convenient instruments for analysis




Recorder's technical specifications	
ECG channels	3
Dimensions (mm)	67*62*19
Weight with battery	70 g
Display	OLED
Memory Card (microSD)	2 Gb
Recording time	24h/48h/72h
Power supply	1 AAA battery
PC connection	Bluetooth

What's included	
ECG recorder	1
ECG cable	1
Charger	1
Rechargeable AAA battery	2
Memory card	2
Bluetooth PC adapter	1
Card reader	1
Case	1


Software CloudHolter

The Holter record can be analyzed using the desktop software or using Cardiomatics AI cloud interpretation service

Desktop software



• Just record, upload and get automatic report



- HR, RR, ST, QT trends
- Automatic QRS recognition
- Automatic ST -segment analysis
- HRV analysis
- Automatic arrhythmia detection
- Pacemaker detection and analysis
- Automatic generation of final report and user templates

Oy Diagnostic Devices Development
Dunckerinkatu 5 AS 1
70500 Kuopio, Finland
www.ddd-medical.com
info@ddd-medical.com