


# V10 TABLET VETERINARY DIAGNOSTIC ULTRASOUND SYSTEM

V10's sleek and stylish curves give emphasis to modern industrial design philosophy. The device has a full touch screen and can be charged quickly. V10 veterinary ultrasound system has a large battery capacity, supporting more than 4 hours of working time, and demonstrates excellent versatility in different fields.




*Charge 2 Hours  
Work 4 Hours*



## SYSTEM CONFIGURATION

Display	8" high resolution LCD Monitor, Full Touch Sreen
Start-up time	≤30 second
Dimension	210mm(L) x 142mm(W) x 16mm(H)
Weight	600g
Connector	1 standard transducer connector
Imaging Modes	B, 2B, B+M, M, 4B,
Port	Blue tooth, Wi-Fi, USB3.0, HDMI, Audio
Battery	Internal Li Ion Battery pack, ≥4 hrs
OB Measurement	Swine
	Bovine
	Equine
	Ovine
	Feline
	Canine
Optional Accessories	Silicon Cover
	Carrying Bag
	Car Charger
	S-Video, VGA, DVI, Type-C
	Transducers
	Sun Shade
	Neck Strap



## TRANSDUCERS



Linear Rectal



Micro-convex



Linear



Convex