



EN-1029

**HD Image**

Flexible Video  
Rhinolaryngoscope



### Multi-function Buttons

Capture, record, freeze, zoom in/out.



### Diagnosis Type

Ultra slim insertion tube for pediatric examination and diagnosis.


### CMOS + LED Technology

CMOS image sensor delivers a high quality image. No fiber: durable and easy maintenance.



### Portable

Compatible with portable monitor HV-3101.



### Application

Rhinolaryngological examination and diagnosis in E.N.T.



Model	Front Tip (mm)	Insertion Tube(mm)	Channel (mm)	Working Length(mm)	Bending (Up/Down)	Field of View
EN-1029	2.5	2.2	/	300	130°/130°	90°