



SEESHEEN
see the endo-sheen

CMOS+LED

FLEXIBLE VIDEO
RHINOLARYNGOSCOPE





USB Output

Compatible to hardware with USB video input. Everywhere E.N.T.

CMOS + LED Technology

CMOS image sensor delivers a high quality image.
No fiber: durable and easy maintenance.



HD Image

High definition image enable more detailed observation of nasal cavity, pharynx and larynx.



Multi-function Buttons

White balance, brightness, freeze

Application

Rhinolaryngological examination and diagnosis in E.N.T.



Slim Insertion Tube
3.4mm rhinolaryngoscope

Model	Insertion Tube(mm)	Channel (mm)	Working Length(mm)	Bending (Up/Down)	Field of View	Depth of Field(mm)
EN-E19H	3.4	/	300	140°/140°	100°	3-100