

CE 0197



**SEESHEEN**  
see the endo-sheen

# CMOS+LED

FLEXIBLE VIDEO  
RHINOLARYNGOSCOPE



EN-1249



### Multi-function Buttons

Capture, record, freeze, zoom in/out.




### Therapy Type

2.0mm channel for interventional treatment.


### CMOS + LED Technology

CMOS image sensor delivers a high quality image.  
No fiber: durable and easy maintenance.



### Portable

Compatible with portable monitor CV-1080/HV-3101.



### Application

Rhinolaryngological examination, diagnosis and therapy in E.N.T.



Model	Insertion Tube(mm)	Channel (mm)	Working Length(mm)	Bending (Up/Down)	Field of View	Depth of Field(mm)
EN-1249	4.8	2.0	365	130°/130°	90°	3-50