

# SKE005

## TRANSFUSION CHAIR



# SKE005

## FEATURES AND HIGHLIGHTS

### Technical parameters

▪ Seat Size (LxWxH)	750×1000×1050 mm
▪ Bed Size (LxWxH)	750×1550×450 mm

### Technical configuration

<input checked="" type="checkbox"/> Anti-slip Feet Cover	4pcs
<input checked="" type="checkbox"/> Wooden Armrest	2pcs
<input checked="" type="checkbox"/> IV Pole Previsions	2pcs
<input checked="" type="checkbox"/> IV Pole	1set
<input checked="" type="checkbox"/> Utility Basket	1pc



#### Design

Scientific design, ergonomics, simple structure, convenient operation, applied to infusion and patient rest. Can pull out to be a bed.



#### Material

Plastic steel frame, strong stability, durable use, beautiful appearance, comfortable texture, convenient cleaning.



#### Back and Seat Plate

The back plate and seat plate are made of high-quality artificial leather and sponge, providing a safe, comfortable and warm infusion experience.




#### Armrest

The handrail is made of imported mountain wood, resin paint, and the shape is generous.



#### Footrest

The pedal board is welded with high-quality steel pipes, which is well-made, sturdy and durable, and has strong rust resistance and fire resistance.



#### Shelf

There is a shelf on the side of the chair, which is convenient for the infusion person to put their belongings. Carbon steel structure, strong bearing capacity.



#### Foot Pad

Beautiful appearance, easy to use and has the functions of anti-slip, anti-collision, anti-shock, temperature resistance, noise reduction, etc.



#### IV Pole

Full stainless steel IV pole, bearing 15 kg.





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

# SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎️ +86 18021231901 · 🕒 24/7

✉️ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue