



## HI-dop

Fast and reliable vascular doppler

- LCD display of pulse frequency and device status
- Available 4, 5 and 8MHz probes
- External volume control
- Connector for earphone
- Automatic "power off" function
- Battery operation with 2×AA size



### Technical Specification

- Pulse rate of 50 ~ 240 beats per minute
- Accuracy  $\pm 2\%$
- Ultrasound probes with 4MHz, 5MHz and 8MHz
- Ultrasound intensity  $< 10\text{mW}/\text{cm}^2$
- Speaker : 1.2W(built-in loud speaker)
- 3.5mm connector for earphone
- Automatic power-off after 5 minutes
- Battery operated by 2 AA batteries (1.5V)
- Battery lasts about 360 minutes
- Including bag, batteries and ultrasound gel

### 3 Probes



3 compatible probes available : 4MHz, 5MHz and 8MHz

### Connector for earphone



3.5mm connector for earphone or external recorder

### Standard Accessories



medical **ECONET**

Mainstrasse 6 c-d, D-45768, Marl, Germany  
 TEL. +49 (0) 2365 924 370  
 FAX. +49 (0) 2365 924 3755  
 Email info@medical-econet.com  
 www.medical-econet.com

CE 0470

distributor

\*All specifications are subject to change without notice.