

INFRARED Thera

A heating lamp mounted on a flexible gooseneck (86 cm long, 34 cm of which are rigid tube). Equipped with a handle, it is easily positioned at the required distance. The model with a mechanical timer switch enables the length of treatment to be set from zero to 60 minutes with an audible signal upon stopping.

Use : Place the infrared lamp at a distance of around 40 cm: the degree of heat must always remain agreeable. Sessions of 10 to 30 minutes are sufficient.

Ideal for physiotherapy, osteopathy and chiropractic use

● Medical device class 1

- Wavelength: 1200 nm
- Electric consumption: 250 Watts
- IR (E27) bulb life span: 5000 h minimum
- Electrical protection: insulation class I
- Standards: IEC 60 601 -1/ -2-41
- Power cable length: 4 meters
- Suggested colors: white
- Head dimensions: Diam 26 cm

IR bulb **87115005**



Thera
Timer switch



- Fixings to choose from p23-26



CODES			
	On/off Switch	Contactless switch	Timer switch
Flexible 86 cm	IRT250860	IRT250860-ISC	IRT250860MI



**HIGH
HEAT OUTPUT**