

Hepta^{LED}



Professional examination lamp with a timeless design that combines power, solidity, and versatility for heavy use.

With a tough flexible neck and an ergonomic handle, it can be positioned as required and fulfils all the needs of medical lighting.

The optional dimmer switch allows the brightness to be adjusted for user comfort and optimal use.

Ideal for all types of use in general medicine, dermatology, gynaecology, urology, podiatry, etc.


● Medical device class 1

Light intensity **980 Lumens***

Colour of temperature **4500°K* (neutral white)**

Colour rendering index **CRI 95/100**

Electric consumption **7 Watts**

Lux*	Beam 10°	Distance	Illuminated diameter
130 000		30 cm	10 cm
60 000		50 cm	12 cm
18 000		100 cm	20 cm

* Values accurate within +/- 10 %



- LED life span: 50 000 h minimum
- Electrical protection: insulation class II
- Standards: IEC 60 601 -1/ -2-41
- Power cable length: 3 meters
- Suggested colors: white and black
- Head dimensions: 11 x 10 cm



Hepta 65 |
Dimmer switch

- Fittings: 1, 2, 3, 4, 6, 7, 8, 9, 10, 11, 12, 13

CODES			
	On/off Switch	Dimmer switch	Contactless switch
Flexible 65 cm	LED07651	LED07651V	LED07651-ISC
Flexible 114 cm	LED071001	LED071001V	LED071001-ISC