

Vera

The key tool for any analysis that requires the highest precision and detail.

A LED powered circular magnifying lamp with a 12 cm diameter 2.25 X (5 dioptries) non-distorting mineral glass lens and a retractable protective cover. It provides white homogeneous glare-free light.

The articulated arm and swivel joint enable rapid and precise positioning of the lamp using only one hand.



Perfect for use in general medicine and dermatology

- LED life span: 50 000 h minimum
- Lighting Power at 10 cm : 2500 Lux
- Electrical protection: insulation class I
- Power cable length: 1.5 meters
- Suggested colors: white
- Dimensions: head diam 12 cm articulated arm 86 cm
- Clamp provided

Light intensity **300 Lumens***

Colour of temperature **6400°K* (white)**

Colour rendering index **CRI 95/100**

Electric consumption **4 Watts**

*Values accurate within +/- 10 %

CODES

Lamp (clamp provide, maximum opening 6 cm)	LED00912
Heavy trolley base high 80 cm	PIR00912
Wall bracket	EMU00912