

A powerful lamp on a very robust flexible arm appreciated by specialists who wish to perform small interventions in their practice, and who need skin and tissue tone rendered accurately. The ergonomic handle enables precise and intuitive positioning.

The optional intensity switch (50/100 %) allows the practitioner to view the skin under different light levels.

Ideal for dermatology, and minor or cosmetic surgery


● Medical device class 1

Light intensity **1680 Lumens\***

Colour of temperature **4500°K\* (neutral white)**

Colour rendering index **CRI 95/100**

Electric consumption **17 Watts**

Lux*	Beam 10°	Distance	Illuminated diameter
210 000		30 cm	10 cm
120 000		50 cm	13 cm
33 000		100 cm	20 cm

\* Values accurate within +/- 10 %

Bella  
contactless switch

- LED life span: 50 000 h minimum
- Electrical protection: insulation class II
- Standards: IEC 60 601 -1/ -2-41
- Power cable length: 3 meters
- Suggested colors: white
- Head dimensions: 16 x 13 cm
- Fittings: 1, 2, 3, 4, 6, 7, 8, 9, 11, 12, 13

CODES		
	On/off Switch	Contactless switch
Flexible 70 cm	LED17650	LED17650-ISC
Flexible 70 cm Double intensity option	LED17650DI	LED17650DI-ISC

