



**PATIENT LIFTER - 140 KG  
COD. SOL2**

Electric patient lifter with rechargeable batteries.

Powder-coated steel tubular frame.

Anti-trauma push handle, rounded edges and corners.

Lifting and descending movements operated by electric motor.

Low battery alarm.

Complete with emergency push-button in order to immediately stop the motion and release the up/down piston stroke.

Adjustable base by manual lever.

On 4 castors diam. 100mm (rear castors are lockable).

Equipped with battery charger and sling holder with balancing bar and 2 hooks.

Complete with standard sling made of non-toxic and hypoallergenic polyester.

The sling mesh is easy to wash (up to 95°) and dry.

Patient's weight: 40-90 kg

Load capacity: 140 kg.

Dim.1190x730x630/1510 h mm

**Position 1**

Rev. 00 del 03/02/2021

TECHNICAL SHEETS - - NR. /