GYNAECOLOGY COUCH COD. ht635

Structure made of epoxy-painted steel, easy to clean and disinfect. Padded mattress is covered with fireproof leatherette, available in a variety of colors.

Independent movements of the ergonomic backrest and legrest by electric motor. Trendelenburg adjustment and height adjustment by electric motor. Height adjustment by electric motor.

The couch is equipped with ergonomic and removable padded armrests and pushing handle behind the backrest. Armrests can be tilted, adjusted horizontally, in height and angle; they are synchronized with backrest and provided with lock and unlock device.

Equipped with adjustable knee crutches and removable stainless steel tray, dim. cm 16x16x10 h.

Base on 4 twin swivel castors, 2 lockable, diam. 100 mm

Technical data:

- backrest dim.: 880x540 mm adjustment 0° 80°
- armrests dim: 150x500 mm
- seat dim.: 380x540 mmm;
- legrest dim.: 380x520 mm adjustment 0° 90°
- Trendelenburg adj.: 0 -21°
- height adj: 530 850 mm
- maximum load capacity: 180 kg
- weight: 72 kg

Dim. in sitting position: 900x900x1350/1670 h mm Dim. in lying position: 1800x900x530/850 h mm

Dim.1800x900x530/850 h mm

Position 1



Upholsery colorsBlackGreenLight blueBlueGreyRedOrangeYellowBeige

Rev. 00 del 02/11/2021

TECHNICAL SHEETS - - NR. /

