



Upholstery colors



Rev. 00 del 03/02/2021

TECHNICAL SHEETS - - NR. /

MULTIPURPOSE COUCH - HYDR. BACKREST AND LEGREST COD. HT631

Structure made of epoxy-painted steel, easy to clean and disinfect.
Padded mattress is covered with fireproof leatherette, available in a variety of colors.

Independent movements of the ergonomic backrest and legrest by gas spring.
Trendelenburg adjustment by gas spring.
Height adjustment by electric motor.

The couch is equipped with ergonomic and removable padded armrests and push handle mounted on backrest. Armrests can be tilted, adjusted horizontally, in height and angle; they are synchronized with backrest and provided with lock and unlock device.

Equipped with pushing handle behind the backrest.

Base on 2 fixed front castors and 2 rear lockable castors, diam. 100 mm

Technical data:

- backrest dim.: 880x540 mm - adjustment 0° - 80°
- seat dim.: 470x540 mm;
- legrest dim.: 380x520 mm - adjustment 0° - 90°
- armrests dim. 150x450 mm
- Trendelenburg adj.: 0 -22°
- height adj: 530 – 850 mm
- maximum load capacity: 180 kg
- weight: 72 kg

Dim. in sitting position: 900x900x1350/1670 h mm

Dim. in lying position: 1800x900x530/850 h mm

Dim. 1800x900x530/850 h mm

Position 1