ENT COUCH - GAS SPRING BACKREST ADJ. COD. HT900A



Height adjustable ENT couch.

Structure made of epoxy-painted steel, easy to clean and disinfect. Padded backrest and seat, upholstered with fireproof leatherette, resistant to disinfectants, blood and urine, available in a variety of colours.

Backrest adjustment by gas spring; fixed legrest. Electric height adjustment from 650 to 980 mm by foot control. Seat rotates 90° clockwise or anticlockwise. Lifting capacity: 180 kg Motor input: 220/230 V - 50 Hz

It is equipped with removable polyurethane armrests

Dim.1800x660x545/920 h mm

Position 1





Rev. 00 del 03/02/2021

TECHNICAL SHEETS - - NR. /

