EXAMINATION TABLE - FIXED H. - 2 SECTIONS COD. HT70

Welded structure made of epoxy-painted tubular steel.

Two-section mattress made of high density polyurethane foam, upholstered with fireproof and waterproof leatherette, available in a variety of colours (standard colour: black)

Leatherette fire resistance class: 1IM.

Manually adjustable backrest.

Paper roll holder available as accessory.

Load capacity: 150kg.

Dim.1800x600x700 h mm

Position 1





Rev. 00 del 03/02/2021

TECHNICAL SHEETS - - NR. /

