



**EXAMINATION TABLE - ELECTRIC ADJ. H. - 2 SECTIONS
COD. HT72**

Structure made of epoxy-painted steel, easy to clean and disinfect.

Two-section mattress made of high density polyurethane foam, upholstered with fireproof and waterproof leatherette, available in a variety of colours. Mattress thickness: 70 mm. Leatherette fire resistance class: 1IM.

Manually adjustable backrest from -14° to $+70^{\circ}$

Height adjustment from 520 to 950 mm by electric actuators 24V controlled by handset.

Base on 2 castors diam. 100 mm and 2 feet.

Load capacity: 200 kg.

Dim.1950x650x520/950 h mm

Position 1

Upholstery colors



Rev. 00 del 03/02/2021

TECHNICAL SHEETS - - NR. /

