



## STRETCHER FOR MRI COD. LO503-120

The transport stretcher is used as a patient transport stretcher in the field of magnetic resonance (MR).

Our stretcher is built in stainless steel suitable for MR systems up to 3 Tesla.

The transport stretcher LO503-120 is supplied with the side rails.

Equipped with four smooth-running wheels ( $\varnothing$  100 mm) (all four can be locked in the direction of travel and rotation), which ensure effortless patient transport.

Standard height 650 mm, on request the transport stretcher is available in heights of 700, 750 and 800 mm at no additional cost (please indicate the desired height when ordering).

Thickness of the covering of the head section and of the lying surface 64 mm.  
Head section continuously adjustable up to + 40 ° by means of a locking bar mechanism.

Dim. 1900x795x650/700/750/800 h mm

### Position 1

Rev. 00 del 03/02/2021

TECHNICAL SHEETS - - NR. /

