



**GYNAECOLOGICAL TABLE - FIXED HEIGHT
COD. HT91**

Welded structure made of epoxy-painted tubular steel.

3-section mattress made of high density polyurethane foam, upholstered with fireproof and waterproof leatherette, available in a variety of colour. (Standard colour: black)

Leatherette fire resistance class: 1IM.

Manually adjustable backrest and legrest.

It is equipped with adjustable knee crutches and removable stainless steel basin.

Load capacity: 150 kg

Dim. 1800x600x700h mm

Position 1

Upholstery colors



Rev. 00 del 03/02/2021

TECHNICAL SHEETS - - NR. /

