



**EL. MULTIPURPOSE COUCH - FIXED H.
COD. HT990A**

Therapy chair designed to assist and facilitate the work of physicians and medical staff. It also increase the patient's comfort considerably and therefore help to relieve tension during treatment.

Frame made of high strength powder-coated steel, white colour RAL 9001. Exchangeable upholstery made of artificial leather, grey colour, resistant to disinfectants, blood and urine; it can be removed easily with a few movements for cleaning or maintenance. Upholstery size: 2100x600x95 h mm

It is equipped with upholstered tilt-up armrests and footrest. Armrests are height and longitudinal adjustable for easy entry and leave; their multilayered soft cushion is completely sealed against liquids and ensures improved comfort. Footrest is manually adjustable of 25 cm.

Independent electric adjustments by 2 motors, controlled by liquid-proof 13 buttons handset, class IP67; it is equipped with automatic buttons for Entry position, Bed position.

Basement on 4 twin castors, diam 75 mm, with total lock.

Technical data:

Seat height (fixed): 650 mm
 Seat tilt adjustment (fixed): 7°
 Armrest size: 570x150x65 h mm
 Backrest adjustment: -7° +80°
 Legrest adjustment: 30°
 Maximum load capacity: 200 kg
 Weight: 63 kg
 Power supply: 100-240V AC, 50/60 Hz - Motors: 24V DC, IPX6

Dim.2150x930x650 h mm

Position 1

Rev. 00 del 03/02/2021

TECHNICAL SHEETS - - NR. /

