



**CARDIOLOGY TABLE
COD. HT100**

Structure made of solid beechwood naturally painted.

Two-section mattress made of high density polyurethane foam, height 80 mm, upholstered with fireproof and waterproof leatherette, standard colour beige.
Leatherette fire resistance class: 1IM.

Manually adjustable backrest.

Load capacity: 150kg
Weight: 30 kg

Dim.1900x700x750 h mm

Position 1

Rev. 00 del 03/02/2021

TECHNICAL SHEETS - - NR. /

