



COLPOSCOPE AC-1000 SERIES




ALLTION®

AC-1320 COLPOSCOPE

(Three-Step, LED, 45° Inclined Viewing)

- Hand polished optics for exceptional clarity
- Individually adjustable binocular optics for true stereoscopic view
- Optical measuring circles for precise measurements
- Clear, bright LED light source that is longer lasting, gives more sure color balance, and is cooler to the touch than the standard halogen light source
- Light intensity >25000 LUX
- Focal length: 300mm
- Magnification: 15x, 7.5x, and 3.75x
- Field of View: 9.5mm, 19mm, 38mm
- Depth of Field: 5mm, 8mm, 32mm
- Smooth adjustable controls for gross and fine focus
- Optic Angle of binocular head is 45°
- Universal Teflon® ball joint never has to be tightened or adjusted and it gives you a fluid 360° motion for gross focus
- Roller base for easy transport
- Space saving small footprint: length is 14mm and width is 12mm
- Clutter reducing internal power cable
- Working height: 33mm to 49mm
- One piece shipping that puts you in operation within minutes of opening the box (no assembly/installation necessary)
- Comes with color instruction booklet, dust cover to limit wear and tear, and two extra fuses



1mm	A	B
3.75X	3.5	17
7.5X	2	8.8
15X	0.9	4.4

Two Circles Above for Measuring Size of Tissue



Items and specification subject to change,
please contact ALLTION for current product information

CE

CE No. 075S019-IT-CE-P02V02
(FDA 510K Number: K061306)
ISO13485 Number: 41532

ALLTION®

Alltion (Guangxi) Instrument Co., Ltd.
Alltion (Wuzhou) Co., Ltd.

Microscopes Are Our Business & Customer Service Is Our Pride

Alltion Building, No.10,3rd Road, Wuzhou Industrial Park, Wuzhou, Guangxi, China. 543000

Tel: +86-774-2836101 (3 lines) Fax: +86-774-2836192

Email: sales@alltion.com or sales@alltion-microscope.com

Web site: www.alltion.com or www.alltion-microscope.com

Printed in March 2016