

CE 0197



SEESTHEEN

# CH Series

FLEXIBLE VIDEO  
CHOLEDOCHOSCOPE



CH1331

CH1349

HD

## Control Button

- ①Zoom
- ②Capture/Recording
- ③Freeze

## CMOS + LED Technology

Fiber-free design makes it more durable.

## Portable

Equipped with image processor HV-3101  
Support external monitor


### CH1349

Insertion tube:  
4.8mm  
Working channel:  
2.0mm



### CH1331

Insertion tube:  
3.1mm  
Working channel:  
1.2mm



## Application

Mainly applied for diagnosis and treatment of gallstones, such as PTCS surgery. The flexible insertion tube can pass a small invasive puncture hole to the common bile duct and hepatobiliary.

Model	Insertion Tube	Channel	Working Length	Bending (Up/Down)	Field of View	Depth of Field	Defintion
CH-1331	3.1mm	1.2mm	450mm	130°/130°	120°	3-50mm	SD
CH-1349	4.8mm	2.0mm	380mm	180°/130°	120°	3-50mm	HD