

COMEN

	Height Adjustable	Fixed Height
Model	BC-107	BC-108
Transparent Baby Basin	✓	✓
Mattress	✓	✓
Brake Casters	✓	✓
Storage Basket	✓	✓
Bed Lift	✓	×
Bed Tilt	✓	✓
Weight	17.5kg	15.0kg
Height Range (Not Including Mattress)	760mm-900mm	860mm
Bed Weight-Bearing	10kg	
Bed Size	860mm*510mm*260mm (L*W*H)	
Applicable Groups	Newborns, infants (0-5 months or babies less than 70cm in length)	



# BC-107/108 Medical Infant Bassinet

Safer sleep. Easier care.



CE



Shenzhen Comen Medical Instruments Co., Ltd.

Add: No.2 of FIYTA Timepiece Building, Nanhuan Avenue, Matian Sub-district,  
Guangming District, Shenzhen, 518106, Guangdong, China.  
Tel:+86-755-2640 8879 400-700-9488 Fax:+86-755-2643 1232  
Website: en.comen.com E-mail: info@szcomen.com

P/N: EN-BC-107/108-4P-20210603-V1.3



## Safe and durable

The structure is stable and reliable. Even it is tilted to the maximum angle of 12, stable structure protects bassinet bed from falling. The lifting/tilting lifespan is above the industry average, which means longer stability and safer operation.

## Skin-friendly and comfortable

Suitable for newborns and babies aged 0-5 months (or less than 70cm in length). The bassinet bed is made of ABS Eco-friendly material, waterproof and breathable, allowing your baby to grow up with peace of mind. Two types of mattress are available to choose. The mattress with waterproof cover is in compliance with biocompatibility standard.



## Strong and trustworthy

To protect baby's safety, we use high-quality carbon steel for the bed frame, which is strong and deformation-resistance.

## Meticulous design

The pressing handle allows the bassinet bed to be tilted at any angle between 0° to 12°, preventing regurgitating milk and reducing mom's worries.



## Clear vision

The transparent bedside walls provide mom a clear, easy view of baby, enabling eye contact. With well design and fine material, it is the best choice for baby care.

## Flexible height

Adjustable height is able to meet different nursing needs.

## Convenience storage

Featuring large space, easy disinfection, the basket allows infant caring within reach.

## Quiet and stable

Medical locking caster wheels enable baby to sleep soundly and sweetly even during transportation.