



AMI | ITALIA
MEDICAL DEVICES

LEARN HOW TO **SAVE A LIFE !**



SMARTLY[®]
SAVER
Mini Trainer

www.amiitalia.com

SMARTCY® SAVER Mini Trainer

CONTROLS AND INDICATORS

Category	Nominal specifications
Buttons	<ul style="list-style-type: none"> ○ ON / OFF (keyboard) ○ ON / OFF (lateral switch) ○ Shock button: to simulate defibrillation shock delivery ○ Adult mode: to select the patient type ○ Paediatric mode: to select the patient type
Visual Indicators	<ul style="list-style-type: none"> ○ "AED ready for use " :1 Green LED ○ "AED Fault /Low battery": 1 Red LED ○ "Place the defibrillation PADS": 2 Red LEDs ○ "Do not touch the patient: 2 Red LEDs ○ "It is possible to touch the patient":1 Green LED ○ Adult patient :1 Green LED ○ Paediatric patient: 1 Green LED ○ ON/OFF button: 2 Green LEDs ○ Shock button: 4 Red LEDs
Sound indicators	Multilingual voices for instructions while using the device Audible warning and danger signals
Speaker	Variable volume through remote control (powered with no. 2 AAA batteries)

PHYSICAL CHARACTERISTICS


Category	Nominal specifications
Dimensions	108x119x42 mm (folded handle) 138x119x42 mm (extended handle)
Weight	280 g (with battery) + 23 g (electrode interconnection cable)

POWERING

Category	Nominal specifications	
Identification code	ST3-0001	
Typology	No. 2 AA batteries	
Voltage / Capacity	3V – 1Ah	
Duration	Stand-by	5 days*
	Continuous operation	6 hours*

*With new and fully charged battery at constant temperature of 20°C and relative humidity without condensation 45%

BOX CONTENT

 <p>MOQ 3 Units</p>	<ul style="list-style-type: none"> ○ 1 Trainer Device ○ 1 Carry Bag ○ 1 Remote Control ○ 1 Pads Cable ○ 1 Pads Set ○ 2 Battery AA (Device) ○ 2 Battery AAA (Remote Control)
--	--