



AMI | ITALIA  
MEDICAL DEVICES

LEARN HOW TO **SAVE A LIFE !**



**SMARTLY**<sup>®</sup>  
**SAVERGEO**  
Mini Trainer

[www.amiitalia.com](http://www.amiitalia.com)

# SMARTLY<sup>®</sup> SAVERGEO Mini Trainer

## CONTROLS AND INDICATORS

Category	Nominal specifications
<b>Buttons</b>	<ul style="list-style-type: none"> <li>○ ON / OFF (keyboard)</li> <li>○ ON / OFF (lateral switch)</li> <li>○ Shock button: to simulate defibrillation shock delivery</li> <li>○ Adult mode: to select the patient type</li> <li>○ Paediatric mode: to select the patient type</li> </ul>
<b>Visual Indicators</b>	<ul style="list-style-type: none"> <li>○ "AED ready for use " :1 Green LED</li> <li>○ "AED Fault /Low battery": 1 Red LED</li> <li>○ "Place the defibrillation PADS": 2 Red LEDs</li> <li>○ "Do not touch the patient: 2 Red LEDs</li> <li>○ "It is possible to touch the patient":1 Green LED</li> <li>○ Adult patient :1 Green LED</li> <li>○ Paediatric patient: 1 Green LED</li> <li>○ ON/OFF button: 2 Green LEDs</li> <li>○ Shock button: 4 Red LEDs</li> <li>○ CPR FEEDBACK 1 Red LED, 1 Orange LED, 1 Yellow LED, 1 Green LED</li> <li>○ CPR bluetooth connection 1 Green LED</li> </ul>
<b>Sound indicators</b>	Multilingual voices for instructions while using the device Audible warning and danger signals
<b>Speaker</b>	Variable volumethrough remote control(powered with no. 2 AAA batteries)

## PHYSICAL CHARACTERISTICS


Category	Nominal specifications
<b>Dimensions</b>	108x119x42 mm (folded handle) 138x119x42 mm (extended handle)
<b>Weight</b>	280 g (with battery) + 23g (electrode interconnection cable)

## POWERING

Category	Nominal specifications	
<b>Identification code</b>	<b>ST7-0001</b>	
<b>Typology</b>	No. 2 AA batteries	
<b>Voltage / Capacity</b>	3V – 1Ah	
<b>Duration</b>	Stand-by	5 days*
	Continuous operation	6 hours*

\*With new and fully charged battery at constant temperature of 20°C and relative humidity without condensation 45%

## BOX CONTENT

 <p>MOQ 3 Units</p>	<ul style="list-style-type: none"> <li>○ 1 Trainer Device</li> <li>○ 1 Carry Bag</li> <li>○ 1 Remote Control</li> <li>○ 1 Pads Cable</li> <li>○ 1 Pads Set</li> <li>○ 2 Battery AA (Device)</li> <li>○ 2 Battery AAA (Remote Control)</li> </ul>
--	--