

Micomme i Series

T3-T5-P1



Your Smart Respiratory Therapy



Sparable design and one button control



CPAP, S, T, S/T, VAT, PC



AST-automatic synchronization technology



COMF-pressure relief technology



VAT-volume assured technology



Preheat and Auto Drying

Micomme i Series



Modes	T3	T5	P1
Mode	CPAP,S,T,S/T	CPAP,S,T,S/T,VAT	CPAP,S,T,S/T,VAT,PC
Pressure range and settings			
Pressure range	4-30cm H2O		
Pressure Accuracy	±2cm H2O		
Max. output pressure	30cm H2O		
Pressure release(COMF)	OFF, 1 - 3 levels		
Ramp time	0min-45min(5-minute increments)		
Humidification Level	1-5 Levels		
Rise Time	1-6Levels		
Air Leak Compensation	60L/min		
Tidal volume	0 - 2000mL		
Maximum flow rate	100L/min@10cmH ₂ O,130L/min@20cmH ₂ O		
Humidification			
Humidification Level	1-5 Levels (95to 167°F/35-75°C)		
Preheat	5 minutes from 23 °C to 45 °C		
Auto Drying	YES		
Data Management			
Data Transmission	USB Key and WiFi		
Dimensions			
Weight	2.1Kg		
Outline Dimensions	30.0cm(L)*21.0cm(W)*13.5cm(H)		
Water Chamber Capacity	320mL±10mL		
Working Voltage	24VDC		