



GIMA

GIMATEMP NO CONTACT INFRARED THERMOMETER IT,DE,PT,SE

Code: 25584
Category: No contact thermometers
Unit of sale: 1 pc.
Minimum order: 1
Type: Medical device
Class: II A
NSIS: 1951667
CND: V03010102
EAN13: 8023279255843



Description: GIMATEMP NON CONTACT THERMOMETER IT, DE, PT, SE

Multifunctional infrared non contact thermometer with 4 applications: body, liquid, food, ambient temperature.

- high accuracy and quick measurement in =2-3 sec
- 3 colour backlight alarm (red, yellow, green) for checking the temperature status
- measuring distance 1-5 cm
- room temperature is always displayed.
- 50 memories

Technical Specifications: Measurement method: Temporal artery scan or forehead spot
Display model - resolution: °C/°F - 0.1°C/0.1°F
Operating mode: Adjusted mode (Body mode)
Direct mode (Surface mode)
Output range:
- Body mode: 34.0~43.0°C/93.2-109.4°F
- Surface mode: 0~100°C/32-212°F
Accuracy: Body mode: 35.0~42.0°C ±0.2°C / 95.0-107.6°F ±0.4°F
±0.3°C/±0.5°F for other range
Surface mode: ±2°C/±3.6°F
Auto power off time: =18 s
Measuring time: 2-3 s
Measuring distance: 1-5 cm (0.4-2 inch)
Batteries (included): 1.5 V (AAA) alkaline batteryX2 (IEC type LR03)
Size - Weight: 122x96x42 mm - 80 g