



SEE CLEARLY TREAT SIMPLY



No capital investment



See and treat in office



No risk of cross-contamination



4.5mm OD, No anesthesia and speculum needed



Compatible with all 5Fr, 7Fr, 9Fr flexible/rigid/semi-rigid instruments




Intuitive orientation in uterus, excellent fluid circulation, and effortless photo/video recordings




MODEL	HTx 40	HTx 60	HTx 60s	HTx 60T
Outer Diameter	4.5mm	6.2mm	6.2mm	6.2mm
Working Channel	2.1mm	3.3mm	3.3mm	3.3mm
Instruments Compatibility	5Fr flexible/semi-rigid	7Fr, 9Fr rigid	MyoSure Lite/Reach/Manual	TruClear mini
Working Length	240mm	240mm	195mm	240mm
Field of View	120°	120°	120°	120°
Direction of View	8°	12°	0°	0°



Tightly sealed without any leakage with tools



Separate in and out water channels for better fluid management



350 degrees knob rotation for complete uterus visibility



Exceptional image quality that fulfills both diagnostic and therapeutic purposes