Relentless Innovation
for your diagnostic confidence

SAMSUNG

V5 ONE

From basics to brilliance



From basics to brilliance

Step into a new era of healthcare with the Samsung V5 ONE ultrasound system, combining outstanding imaging performance and advanced AI diagnostic features. The V5 ONE supports a broad range of medical specialties with superior image clarity. Intelligent AI functionalities streamline diagnostics, enabling quicker, accurate assessments for general imaging, women's health, cardiovascular, and musculoskeletal examinations.

Despite its slim and compact design, the V5 ONE does not compromise on power, offering exceptional mobility and usability. Its intuitive controls and automated features enhance workflow efficiency, simplifying daily operations and setting new standards in ultrasound diagnostics. The Samsung V5 ONE is designed to meet the demanding needs of modern healthcare professionals, making it an indispensable tool in any medical setting.

* V5 ONE ultrasound system supports Windows 11 operating system and complies with the latest security standards.



Learn more


Striking images for extra confidence




Liver Hemangioma




Renal in MV-Flow™[1]



Trivial TR in modified view



Suprapatellar tendon with HQ-Vision™



Endometrial Polyp



Embryo in RealisticVue™[1]


1. Optional feature, additional purchase required.

Elevate diagnostic precision in general imaging

Transform your routine ultrasound assessments with the V5 ONE ultrasound system, engineered for comprehensive support across general imaging, cardiovascular, and musculoskeletal diagnostics. Harness the power of our advanced automation technologies that optimize your diagnostic process, enabling you to deliver precise, dependable results with ease.

Display and quantify tissue stiffness in a non-invasive method

S-Shearwave Imaging™[1] allows for non-invasive assessment of the stiffness of tissue/lesions in various applications such as breast and liver. The color-coded elastogram, quantitative measurements, dual or single display option, and user-selectable ROI functions are especially useful for the accurate diagnosis of breast and liver diseases.



Quantitative measurement of liver fat with ultrasound signal

TAI™[1] (Tissue Attenuation Imaging) provides quantitative tissue attenuation measurement to assess steatotic liver changes.

TSI™[1] (Tissue Scatter distribution Imaging) provides quantitative tissue scatter distribution measurement to assess steatotic liver changes.




Measure the size of the organ of interest automatically and simplify your workflows

* Target : Kidney, Spleen, Bladder, Prostate, Bowel


Kidney in AbdomenAssist™[1] - Reduced keystrokes by about **75%** compared to manual input



BladderAssist™[1] - Reduced keystrokes by about **89%** compared to manual input



ProstateAssist™[1] - Reduced keystrokes by about **89%** compared to manual input




1. Optional feature, additional purchase required.

Easy calculation of the strain ratio between two ROIs

E-Strain™[1, 2] is designed to enable quick and easy calculation of the strain ratio between two regions of interest for day-to-day practice. Simply by setting the two targets, you can receive accurate, consistent results and make informed decisions in many types of diagnostic procedures.

Quantify wall motion of the left ventricle


Strain+[1] is a quantitative tool for measuring global and segmental wall motion of the left ventricle (LV). Three standard LV views and a Bull's Eye are displayed in a quad screen for easy assessment of the LV function.



Analyze selected breast lesions and report breast assessment


S-Detect™[1, 3] *for breast*, analyzes selected lesions in the breast ultrasound study and shows the analysis data, applies BI-RADS ATLAS* (Breast Imaging-Reporting and Data System, Atlas) to provide standardized reporting; and helps diagnosis with the streamlined workflow.

* It is a registered trademark of ACR and all rights reserved by ACR.



Detect and track nerves automatically with AI

NerveTrack™[1], a feature based on Deep Learning technology, detects and provides information of the location of the nerve area in real-time during ultrasound scanning.




An automated reporting tool for heart diagnosis

HeartAssist™[1], a feature based on Deep Learning technology, provides automatic classification of ultrasound image into measurement views required for heart diagnosis and provides measurement results.

Reduced keystrokes by about 90% compared to manual input

* 2D - Simpson 4ch view




1. Optional feature, additional purchase required.

Unlock comprehensive women's health solutions

Enhance your daily women's health assessments with the V5 ONE ultrasound system, tailored for unmatched support in obstetrics and gynecology. This system integrates our most advanced automation technologies to streamline your diagnostic process in women's health, ensuring precise and reliable results effortlessly.


Assess the risk of infertility using volume data

5D Follicle™[1] identifies and measures multiple ovarian follicles in one scan for rapid assessment of follicular size and status during controlled ovarian stimulation.




Automated support for freezing, annotation, and measurement of fetal structures in real-time

Live ViewAssist™[1], a feature based on Deep Learning technology, automatically classifies ultrasound images in real-time and provides annotation of structures and measurement results.




Support clinical standards with AI

QualityCheck is based on AI technology to assess whether the views acquired during prenatal ultrasound examinations meet clinical standards. It provides Quality Criteria to assist users in detecting with appropriate criteria and maintaining image quality consistency.



Measure the size and shape of the uterus with AI technology


UterineAssist™[1], based on Deep Learning technology, automatically measures the size and shape of the uterus, assisting in detecting signs of uterine-related abnormalities, as well as reducing scan time.



1. Optional feature, additional purchase required.


Enhance daily efficiency

The V5 ONE ultrasound system enhances efficiency with advanced workflow solutions, remote accessibility, a larger screen view, and a compact, robust structure with convenient battery. These elements make it ideal for healthcare professionals seeking to improve productivity and patient care.



See images in expanded view

The ultrasound examination can be performed while viewing the images and cines that are expanded at various ratios according to the user preference.




Real-time image sharing solution

SonoSync™[1, 4] is available in PC and smartphone, etc. as a real-time image share solution that allows communication for care guide and training between doctors and sonographers. In addition, voice chatting, text chatting, video conference, and real-time marking functions are provided for better communication; and the MultiVue function is included that allows monitoring multiple ultrasound images on a single screen.


Continue working even when AC power is temporarily unavailable

BatteryAssist™[1] provides battery power to the system, enabling users to perform scans when AC power is temporarily unavailable. It also allows to transport the ultrasound system to another location and start to scan right away.



Build predefined protocols to ensure every step is followed every time

EzExam+™ enables you to build or use a predefined protocol, and assign protocols for examinations that are regularly performed in the hospital in order to reduce the number of steps that you have to go through.





1. Optional feature, additional purchase required.

Comprehensive selection of transducers

Curved array transducers



CA1-7SD *

Abdomen, Obstetrics, Gynecology, Pediatric, Musculoskeletal, Vascular, Urology, Thoracic



CA2-8AD

Abdomen, Obstetrics, Gynecology, Pediatric, Musculoskeletal, Vascular, Urology



CA4-10M *

Abdomen, Pediatric, Vascular

Endocavity transducers



EA2-11ARE *

Obstetrics, Gynecology, Urology



EA2-11AVE *

Obstetrics, Gynecology, Urology



miniER7 *

Obstetrics, Gynecology, Urology

Volume transducers



CV1-8AE

Abdomen, Obstetrics, Gynecology, Urology



EV2-10A *

Obstetrics, Gynecology, Urology

Linear array transducers



LA2-9S *

Abdomen, Pediatric, Musculoskeletal, Vascular, Small parts



LA2-16S

Abdomen, Pediatric, Musculoskeletal, Vascular, Small parts



LA2-16Sc

Abdomen, Pediatric, Musculoskeletal, Vascular, Small parts



LA3-14AD

Abdomen, Pediatric, Musculoskeletal, Vascular, Small parts



L3-22

Musculoskeletal, Pediatric, Vascular, Small parts, Dermatology



LA3-22AI

Musculoskeletal, Intraoperative, Dermatology

Phased array transducers



PA1-5AE *

Cardiac, Vascular, Abdomen, Pediatric, TCD



PA2-9S

Cardiac, Vascular, Abdomen, Pediatric, TCD

CW transducers



DP2B

Cardiac, Vascular, TCD



CW6.0

Cardiac, Vascular, TCD

TEE transducers



MMPT3-7

Cardiac



TA2-9

Cardiac

Bi-plane transducer



BCC2-13

Urology

* Ergonomic transducers
The new endocavity transducer supports natural grip by moving the max-width point to a more forward position and also increasing the length of the grip to allow balanced weight distribution.

Cleaning and disinfection guide

1. Optional feature which may require additional purchase.
2. Strain value for ElastoScan+™ is not applicable in the United States and Canada.
3. Recommendations about whether results are benign or malignant in S-Detect™ are not applicable in the United States.
4. SonoSync™ is a function for image sharing, not for diagnosis.

* This product, features, options, and transducers may not be commercially available in some countries.
* Sales and Shipments are effective only after the approval by the regulatory affairs.
Please contact your local sales representative for further details.
* This product is a medical device, please read the user manual carefully before use.
* S-Vue Transducer™ is the name of Samsung's advanced transducer technology.
* V5 ONE ultrasound system supports Windows 11 operating system and complies with the latest security standards.
* Keystroke reduction percentages are based on internal data comparing AI-assisted features to manual measurement methods.

SAMSUNG MEDISON CO., LTD.

© 2026 Samsung Medison All Rights Reserved.
Samsung Medison reserves the right to modify the design, packaging, specifications, and features shown herein, without prior notice or obligation.

CE 0123

CT-V5-v2.01-ALL-DESIGNSTICK-261002-EN